


# Your Siemens digital hearing instrument

All digital hearing instrument are different, some have volume controls and buttons, others do not. Your audiologist will inform you about the controls on your specific hearing instrument(s).

## Features:

-  Ear hook
  -  Microphones
  -  Volume control
  -  Program/on and off button
  -  Battery compartment
  -  Colour plate
- Red = Right hearing instrument  
Blue = Left hearing instrument



### Siemens Hearing Instruments Ltd.

Alexandra House, Newton Road, Manor Royal, Crawley, West Sussex RH10 9TT

Tel: 01293 423700 Fax: 01293 403080


Email: [siemenshearing.shi.uk@siemens.com](mailto:siemenshearing.shi.uk@siemens.com)

[www.siemens.co.uk/hearing](http://www.siemens.co.uk/hearing)

**SIEMENS**



## Changing the battery

1. Open the battery compartment.



2. Remove coloured sticker from new battery, place battery in drawer with + symbol on battery + side up. (You should only remove the protective sticker when you actually use the battery).


3. Close battery door.



## On and off


To turn the hearing instrument on, close battery compartment.

**on**



To turn the hearing instrument off, partially open the battery compartment.

**off**



## The program button

The button is used to change the hearing instrument program.



Your audiologist will complete the details:

- |  |        |                          |          |                          |
|--|--------|--------------------------|----------|--------------------------|
| 1. Everyday listening = You will hear 1 beep | Active | <input type="checkbox"/> | Inactive | <input type="checkbox"/> |
| 2. .... = You will hear 2 beeps              | Active | <input type="checkbox"/> | Inactive | <input type="checkbox"/> |
| 3. .... = You will hear 3 beeps              | Active | <input type="checkbox"/> | Inactive | <input type="checkbox"/> |

---

## Using the button to switch on and off

Can I use the button to turn  
the hearing instrument on and off?    Yes     No

- To turn the hearing instrument off:  
press and hold the button for 3 seconds.
- To turn the hearing instrument on:  
press and hold the button for 3 seconds.

**Note:** When the hearing instrument is turned off using the button, wait 10 seconds before trying to turn the instrument on again.


## The volume control

Manual control can be fitted on some instruments?

Can I use the manual volume control? Yes  No

To adjust the volume:

- Roll the wheel up to make it louder.
- Roll the wheel down to make it quieter.



## Using the telephone with AutoPhone™

The AutoPhone function is:

Active  Inactive

- Hold the telephone next to your hearing instrument at the top of your ear. The hearing instrument will automatically change to AutoPhone and then return to normal when the telephone is moved away.



## Cleaning your hearing instrument

Your audiologist will give you full explanation of how to clean your hearing instrument, however, you may use a damp cloth to wipe over the hearing instrument and ear mould.

